


# Blackridge GPS Watch

## WANDERING PREVENTION



# Safety for Alzheimer's, Dementia & Autism

## Alerting and GPS locating system for those who wander

Designed specifically for the "at-risk" market, in particular those that are prone to wandering (i.e.: Alzheimer's, Dementia & Autism), the Blackridge GPS Watch provides caregivers and loved ones a greater degree of safety and peace of mind, day and night.

Introduced to the market in late 2016, the Blackridge GPS Watch is not only the smallest, lightest and most aesthetically pleasing GPS Watch technology on the market, but it is fully loaded with features which allow for greater freedom and mobility to those who may be at risk to wandering.

The ability to have a two-way voice conversation directly to the watch allows for immediate communication, while at the same time providing a highly accurate GPS locate every minute. Viewable on Google maps on any iPhone or Android smartphone, care providers will know the exact location of the GPS Watch at all times.

Worried about the watch being removed? The Blackridge GPS Watch allows for the watch strap to be locked so that a special tool is required for removal. Also, the watch is equipped with an alert removal notification, which will allow care providers to be aware when the watch has been removed, if the lock strap option is not chosen.

Other common challenges, such as battery life and waterproof have also been addressed by the Blackridge GPS Watch. Go up to 4 days in between charges and apply our waterproof shield to bring the watch up to an IP66 rating.

Finally technology has caught up to wandering prevention requirements and is now available at a low cost from BLACKRIDGE SOLUTIONS.

**"When safety and peace of mind is important"**



## Two-way Voice Communication

The watch has the ability to make and receive cellular based voice calls

## Lockable Watch Strap

Requires special tool to unlock. Ensures the watch can not be removed by the individual, if preferred

## Long Battery Life

Go up to (4) days in-between charges. Battery level is displayed on the front of the watch

## Waterproof

A waterproof shield is available for the watch

## Detailed GPS Tracking (Indoors/Outdoors)

New GPS locates are provided every minute and are viewable on the smartphone app

## Panic Button (SOS)

The watch has a Panic button on the side to allow for an immediate call out to up to (16) pre-defined telephone numbers, along with providing a current GPS locate

## Geofencing

Define one or multiple areas on a map and be made aware automatically when the watch enters or exits the area

## Step Counter

The front display of the watch shows how many steps the user has taken during the day


## Alarm Clock

The watch will provide an audible alarm at specific dates and times, alerting the user of a specific task or event that has been pre-defined

# Blackridge GPS Watch

## WANDERING PREVENTION

TWO WAY VOICE



## FREQUENTLY ASKED QUESTIONS (FAQ)

**Is the cost of the GPS Watch and Service Plan covered by Medicaid or Private Insurance?**

ANSWER: Yes, the GPS Watch and first year Service Plan are covered – USA Only (Inquire for further details)

**How heavy is the GPS Watch and what wrist size will it fit on?**

ANSWER: The GPS Watch is very light ... only 2.2 oz and can fit a wrist between 5 ½" - 8 ½"

**Besides the purchase price for the GPS Watch (\$129.99 USD), are there any other costs involved?**

ANSWER: Yes, a Monthly Service Plan of \$29.95 is required for the GPS data & voice and a one-time \$50.00 Activation Support Fee

**Does the Monthly Service Plan require a minimum commitment period?**

ANSWER: No, the Service Plan is billed monthly to a credit card and can be cancelled at any time at no charge.

**Is the GPS Watch waterproof?**

ANSWER: Yes, with the purchase of our Enhanced Waterproof Shield, the watch is then rated at IP-66 (Heavy spray, no submersion)

**Does the GPS Watch have an option to lock the watch strap so that it can't be removed?**

ANSWER: Yes, with the purchase of our Enhanced Locking Strap, a special tool is required to unlock the watch strap

**How frequently will the GPS Watch provide a GPS locate?**

ANSWER: The GPS locates can be set to automatically report every minute or on-demand

**Does the GPS Watch have an SOS button for emergency situations?**

ANSWER: Yes, by pressing on the SOS button, the GPS Watch will automatically begin calling pre-defined phone numbers

**Can I call the GPS Watch and listen unannounced to the environment where the watch is?**

ANSWER: Yes, the GPS Watch provides a "listen in" feature. This feature can be disabled, if required, for privacy purposes.

**Can I see a historical trail of GPS locates showing where the GPS Watch has been?**

ANSWER: Yes, the smartphone app allows you to view a historical breadcrumb trail, selectable by date and time

**How many care providers can access the secure smartphone app?**

ANSWER: You can assign up to (16) care providers per GPS Watch

**Can the GPS Watch provide entry and exit alert notifications for pre-defined locations (geofences)?**

ANSWER: Yes, you can be notified when the GPS Watch enters and/or exits up to (8) pre-defined locations

**Can the GPS Watch alert me if the user has taken it off?**

ANSWER: Yes, the GPS Watch has a light sensor on the back of it to allow for removal notification

**Can the GPS Watch alert me on pre-defined specific dates and times (ie: medications, appointments, etc.)**

ANSWER: Yes, you can set up dates and times for the GPS Watch to sound a notification reminder

**How will I know when the GPS Watch requires recharging and how frequently is that required?**

ANSWER: The GPS Watch provides a battery level on its display, and it is recommended to be recharged every (4) days

**Can the GPS Watch be powered off by the individual wearing it?**

ANSWER: No, the GPS Watch can only be powered off through the smartphone app

**Can the GPS Watch provide details as to the activity level of the user?**

ANSWER: Yes, the GPS Watch provides a step counter indicator on the front display



**BLACKRIDGE**  
SOLUTIONS